



» **MC Plus**
For R&D

» **CellT Plus**
For cGMP

ADAM™ MC Plus & CellT Plus

Most Accurate Fluorescence Cell Counter


| | |
|---------------|---------------------------------------|
| PBMCs | 10 sec per test |
| Stem cells | 15 µL sample loading volume |
| Primary cells | 3.2 µL measuring volume |
| Cell lines | 13 captured images per channel |

ADAM™ MC Plus & CellT Plus are


the highly accurate automated fluorescence cell counter equipped with bright field and two fluorescent channels (AO/DAPI).

ADAM™ MC Plus is used for R&D, ADAM™ CellT Plus is available in cGMP production environment.

4 tests on 1 slide



Bright field &
Dual fluorescence (AO/DAPI)



▲ AO (Live)


▲ DAPI (Dead)

ADAM™ MC Plus & CellT Plus


ADAM™ MC Plus & CellT Plus measures the number of total cells, viable cells, non-viable cells and shows viability results. In addition, they analyze the cell size and cell aggregation ratio as well.

Applicable to a various cell lines


It is possible to use ADAM™ MC Plus & CellT Plus depending on the cell types (PBMCs, etc.) that need to be monitored during the manufacturing of cell therapy products.




▲ Cell Size




▲ Cell aggregation ratio (Aggregation 36.43%)



PBMCs



Stem cells



Primary cells



Cell lines

PBMCs
Clinical immunology

Stem cells
Regenerative medicine


Primary cells/Cell lines
Basic research

Splenocytes
Vaccine development

Various applications

Large measuring volume

13 images
3.2 µL measuring volume




Large measurement volume obtained by detecting multiple images of the samples through moving stage provides more accurate results compared to other manufacturers' products.


- Cell counting accuracy depends on the counting volume (measurement volume).
- The obtained multi-images are processed by image analysis software integrated inside the system.

Accuracy & Repeatability

Correlation of PBMCs total counting between flow cytometry, manual count and ADAM™ MC Plus & CellT Plus.




Data were compared between flow cytometer and ADAM™ MC Plus & CellT Plus results in 20 different PBMCs concentration samples.



Data were compared between manual count and ADAM™ MC Plus results in 20 different PBMCs concentration samples.

Linearity & Reproducibility



A high-concentration (1.8E+07cells/mL) of PBMCs was diluted and a dilution series was counted by ADAM™ MC Plus & CellT Plus. It shows excellent dilution linearity and reproducibility.

| | High | Medium | Low |
|--------|----------|----------|----------|
| Mean | 1.8.E+07 | 9.0.E+06 | 1.6.E+06 |
| SD | 2.7.E+05 | 2.3.E+05 | 3.8.E+04 |
| CV (%) | 1.5 % | 2.5 % | 2.3 % |

Sample with low, medium and high concentration of PBMCs were counted with three ADAM™ MC Plus & CellT Plus.

Accurate result

Specifications I



ADAM™ MC Plus & CellT Plus

| Hardware | |
|------------------------|--|
| Measuring range | 5×10 ⁴ ~ 2×10 ⁷ cells / mL |
| Optimal range | 4×10 ⁵ ~ 4×10 ⁶ cells / mL |
| Analysis time | fast mode 10 sec / test real cell size mode 30 sec / test |
| Measuring volume | 3.2 µL |
| Focus | Auto-focusing |
| Objective lens | 4 X |
| Weight | 7.0 kg |
| Dimensions (W × D × H) | 277 x 276 x 270 mm |

AccuPlus Slide & Reagent

| Performance | |
|-----------------------|---|
| Staining method | Acridine orange (AO) & 4',6-diamidino-2-phenylindole (DAPI) stain |
| Sample loading volume | 15 µL/test |

Ordering Information

| Cat. No. | Description | Contents |
|-----------------|-------------------------------------|---|
| ADAM-MC Plus | Fluorescence cell analyzer | <ul style="list-style-type: none"> Main device User manual |
| ADAM-CellT Plus | Fluorescence cell analyzer for cGMP | <ul style="list-style-type: none"> Main device User manual 21 CFR PART 11 requirement support appendix |
| APAD-400 | Cell viability reagent | <ul style="list-style-type: none"> 20 mL x 1 bottle |
| AP4S-100 | AccuPlus Slide 4ch. | <ul style="list-style-type: none"> 4ch. Slide 100 ea |
| APB-001 | Test Beads | <ul style="list-style-type: none"> 1 mL x 1 tube / pack |
| QCS-002 | QC Slide | <ul style="list-style-type: none"> 1 QC Slide / case User manual |

FOR RESEARCH USE ONLY. This product is not approved for diagnostic or therapeutic use.



website | www.nanoentek.com
e-mail | sales@nanoentek.com

NESCT-AMCCTP-001E (V.0.1)

NanoEntek, Inc.

Head Office
12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea
Tel +82-2-6220-7940 / Fax +82-2-6220-7999

NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA
Tel +1-781-472-2558 / Fax +1-781-790-5649

NanoEntek Europe | mts med-tech supplies GmbH

Lochhamerstr. 4a, 82152 Martinsried, Germany
Tel +49-89-21-55-38-43 / Fax +49-89-99-95-46-60