



Disposable chamber

for sperm counting

DCS-S01

NanoEnTek

NanoEnTek Inc.

Website : www.nanoentek.com
E-mail : sales@nanoentek.com

Manufactured by
NanoEnTek, Inc.
851-14, Seohae-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do, 18531, Korea
Tel : +82-2-6220-7940 / Fax : +82-2-6220-7999

US Corporation
NanoEnTek America, Inc.
240 Bear Hill Road, Suite 101, Waltham, MA 02451, USA
Tel : +1-781-472-2558 / Fax : +1-781-790-5649

EC REP

MT Promedt Consulting GmbH
Altenhofstrasse, 80, 66386 St. Ingbert, Germany

All the materials in this manual are protected by Korean and international copyright laws.
This manual may not be reproduced, translated, published or distributed by any means without the permission of the copyright owner.

The information in this manual is described as correctly as possible, but may be changed without prior consent or notification.

Documentation : NESMU-DCS-001E (V.0.1)

Revision history: V.0.0 Date: JUL 2014
V.0.1 Date: MAR 2019

Copyright © 2014, by NanoEnTek Inc. All rights reserved. Published in Korea.

Unpacking

When you receive the S-Chip (DCS-S01) for the first time, you will find the following components in your package.

Disposable chamber for sperm counting
Instruction manual

Safety Precautions

For analyzing hazardous or potential infectious materials:

- Take necessary precautions
- Handle with care
- Dispose in an appropriate way

Long exposure to solvents will cause the slide to warp.
Xylene and toluene based mounting media should be avoided.
Glycerol, gelatin, and other aqueous-based media are recommended.

Safety Symbols

The safety symbols on the S-chip (DCS-S01) are intended to inform you of potential danger or a particular caution. Before using, please read and follow the guide for the symbols and their meanings.

LOT Batch code (Lot Number)

Use by

Do not reuse

Manufacturer

Consult instructions for use

NOTE : The S-chip (DCS-S01) is for **single use only**. **Do not reuse**. It should be used immediately after unsealing. The warranty on the S-chip included in the conditions of supply is valid for 24months from the date of manufacture. The **expiration date** is indicated on the front side of outer box.

Counting with S-chip

A. Method

- 1) Mix the specimen thoroughly.
- 2) Load the sample (2 μ L) to the one chamber of S-chip using proper pipette and tip.
- 3) Count the sperm in ten squares under the microscope.

B. Counting

The S-chip consists of 100 small squares.
Each small square measure 0.1 x 0.1 mm, with a chamber depth of 0.02 mm.
Each small square has a total volume of 2.0×10^{-4} mm³.
Count all the sperm in the 10 small squares. The large square has a surface area of 1.0 mm² and chamber depth of 0.02 mm. The large square has 0.02 mm³. As there are 1000 mm³ per mL, the total volume of 10 squares is 0.00002 mL, so that equates to 1/0.00002 mL = 50,000


Sperm per mL = Number of cells in 10 small squares x dilution factor x 50,000

Introduction

The S-chip(DCS-S01) is a disposable chamber for sperm counting. It consists of surface-patterned two enclosed chambers with two ports for sample injection.


Each chamber has 100 small squares, each measuring 0.1 x 0.1 x 0.02 mm = 0.0002mm³. Each square has a total volume of 0.02 μ L (mm³).

Figure 1. Grid pattern



Sample injection area

Figure 2. DCS-S01



NESMU-DCS-001E (V.0.1)



S-Chip

INSTRUCTIONS

Disposable chamber

for sperm counting

DCS-S01

Nano**EnTek**

NanoEnTek Inc.

Website : www.nanoentek.com

E-mail : sales@nanoentek.com



Manufactured by

NanoEnTek, Inc.

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do, 18531, Korea

Tel : +82-2-6220-7940 / Fax : +82-2-6220-7999

US Corporation

NanoEnTek America, Inc.

240 Bear Hill Road, Suite 101, Waltham, MA 02451, USA

Tel: +1-781-472-2558 / Fax: +1-781-790-5649



MT Promedt Consulting GmbH

Altenhofstrasse, 80, 66386 St. Ingbert, Germany

All the materials in this manual are protected by Korean and international copyright laws.

This manual may not be reproduced, translated, published or distributed by any means without the permission of the copyright owner.

The information in this manual is described as correctly as possible, but may be changed without prior consent or notification.

Documentation : NESMU-DCS-001E (V.0.1)

Revision history: V.0.0 Date: JUL 2014

V.0.1 Date: MAR 2019

Unpacking

When you receive the S-Chip (DCS-S01) for the first time, you will find the following components in your package.

Disposable chamber for sperm counting
Instruction manual

Safety Precautions

For analyzing hazardous or potential infectious materials:

- Take necessary precautions
- Handle with care
- Dispose in an appropriate way

Long exposure to solvents will cause the slide to warp.
Xylene and toluene based mounting media should be avoided.
Glycerol, gelatin, and other aqueous-based media are recommended.

Safety Symbols

The safety symbols on the S-chip (DCS-S01) are intended to inform you of potential danger or a particular caution. Before using, please read and follow the guide for the symbols and their meanings.



Batch code (Lot Number)



Use by



Do not reuse



Manufacturer



Consult instructions for use



NOTE : The S-chip (DCS-S01) is for **single use** only. **Do not reuse.**
It should be used immediately after unsealing. The warranty on the S-chip included in the conditions of supply is valid for 24months from the date of manufacture. The **expiration date** is indicated on the front side of outer box.

Introduction

The S-chip(DCS-S01) is a disposable chamber for sperm counting. It consists of surface-patterned two enclosed chambers with two ports for sample injection.

Each chamber has 100 small squares, each measuring $0.1 \times 0.1 \times 0.02 \text{ mm} = 0.0002 \text{ mm}^3$. Each square has a total volume of $0.02 \text{ }\mu\text{L}$ (mm^3).

Figure 1. Grid pattern



Counting with S-chip

A. Method

- 1) Mix the specimen thoroughly.
- 2) Load the sample (2 μL) to the one chamber of S-chip using proper pipette and tip.
- 3) Count the sperm in ten squares under the microscope.

B. Counting

The S-chip consists of 100 small squares.

Each small square measure $0.1 \times 0.1 \text{ mm}$, with a chamber depth of 0.02 mm .

Each small square has a total volume of $2.0 \times 10^{-4} \text{ mm}^3$.

Count all the sperm in the 10 small squares. The large square has a surface area of 1.0 mm^2 and chamber depth of 0.02 mm . The large square has 0.02 mm^3 . As there are 1000 mm^3 per mL, the total volume of 10 squares is 0.00002 mL , so that equates to $1/0.00002 \text{ mL} = 50,000$

$$\text{Sperm per mL} = \text{Number of cells in 10 small squares} \times \text{dilution factor} \times 50,000$$