


All-in-One

residual WBC & RBC counter



ADAM rXBC Series

Blood Cell Counter

ADAM rXBC Series, all-in-one solution for counting residual WBC and RBC!

ADAM™ rXBC Series automatically counts the number of residual leukocytes (rWBCs) and erythrocytes (rRBCs) in blood products using just one device.

ADAM™ rXBC is a stand-alone device that allows single-sample counting. Once the slide with the loaded sample is inserted, all procedures are automated. The device takes 50 seconds to count one sample.

ADAM™ rXBC HT gives you walkaway convenience and handles high throughput requirements to improve further productivity. It takes 50 minutes to count 50 samples.

Automated cell counting eliminates user bias or subjective interpretation that can be found when counting residual leukocytes or erythrocytes using other methods.

ADAM™ rXBC Series can count



Blood Components Analysis

ADAM™ rXBC Series provides analysis of rWBC and rRBC in blood samples, offering higher accuracy for samples with low WBC and RBC counts.



- rWBC residual White Blood Cells



- rWBC residual White Blood Cells
- rRBC residual Red Blood Cells



- rWBC residual White Blood Cells
- rRBC residual Red Blood Cells

Features



Fast
measurement

< 1 min/test



Easy
to use

Simple and
easy procedure





Accurate
and reliable results

Outstanding linearity
& repeatability.



Performance Test - Linearity

Low and high rRBC concentrations in plasma and platelet samples were subsequently counted by ADAM™ rXBC series, which showed outstanding linearity results.

• Plasma Products



• Platelet Products



Precision Test


Plasma samples were prepared at eight different targeted concentrations and platelet samples at five different targeted concentrations. For each concentration, rRBCs were counted 20 times to assess the repeatability of the ADAM™ rXBC series.

rRBCs in Plasma (n=20)			
Target value (cells/μL)	Mean (cells/μL)	STDEV	CV (%)
0 - 1	0.74	-	-
1 - 10	2.51	0.28	11.06
10 - 30	24.63	1.49	6.04
40 - 60	55.18	1.77	3.21
90 - 200	112.08	2.08	1.86
300 - 600	596.34	6.51	1.09
700 - 1000	1192.62	16.49	1.38
4000 - 6000	4147.49	163.96	3.95

rRBCs in Platelet (n=20)			
Target value (cells/μL)	Mean (cells/μL)	STDEV	CV (%)
0 - 1	0.76	-	-
1 - 10	2.60	0.23	8.70
10 - 30	21.37	0.71	3.31
40 - 60	55.27	1.35	2.44
90 - 200	109.52	2.09	1.91

Performance - Comparison to flow cytometry










Refer to the following comparison of residual white blood cell enumeration methods between flow cytometry (Leucocount) and ADAM™ rXBC series in using different amounts of white blood cells aliquots. This comparison was performed using both RBC and platelet samples.



Specifications

Item	ADAM™ rXBC	ADAM™ rXBC HT
Specimen	rWBC: RBC, Platelet concentrates, Plasma rRBC: Platelet concentrates, Plasma	
Measuring range	1 - 100 WBC cells/ μ L 1 - 6,000 RBC cells/ μ L	
Sample volume	100 μ L	
Analysis time	< 1 min/test	
Power	DC 12V, 5A	AC 100-230VAC, 50/60 Hz
Weight	7 kg	120 kg
Dimension	277(W) x 276(D) x 270(H) mm	1020(W) x 730(D) x 670(H) mm

Ordering information

Cat. No.	Product	Contents
ADAM-rXBC	ADAM™ rXBC, Blood cell counter	 Main instrument  Desktop mini PC set
ADAM-rXBC HT	ADAM™ rXBC HT, High-throughput blood cell counter	 Main instrument  Desktop PC set
ADWK-050	ADAM rWBC Kit residual WBC counting kit (for ADAM™ rXBC & HT, 50 tests)	 r-Slide, 50 slides  r-Solution 25 mL  Standard bead solution, 7 ml
ADRK-050	ADAM rRBC kit residual RBC counting kit (for ADAM™ rXBC & HT, 50 tests)	 rRBC Slide, 50 slides  Standard bead solution, 7 ml



website www.nanoentek.com
e-mail ivdst@nanoentek.com

FOR RESEARCH USE ONLY.
This product is not approved for diagnostic or therapeutic use.