
FREND BNP

REF FRBN 025

- English
- German
- Italian
- Spanish
- French
- Portuguese

FREND BNP

Quantitative Assay for B-type Natriuretic Peptide

FRBN 025

For *in vitro* diagnostic use only

1. Intended use

The FREND BNP test is a fluorescent nanoparticle immunoassay designed for *in vitro* quantitative measurement of B-type natriuretic peptide (BNP) in human plasma samples using EDTA as the anticoagulant. The FREND BNP microfluidic flow cartridge is designed for use in the FREND™ System to aid in the diagnosis and assessment of the severity of congestive heart failure. FREND BNP is for professional use only.

2. Summary and explanation of test

Almost six million people in the US are diagnosed with congestive heart failure (CHF), with more than 500,000 new cases diagnosed each year.^[1,2] A chronic, progressive condition where the heart does not meet the blood supply demands of the body, CHF can occur at all ages but becomes more prevalent in patients over 65 years of age. Common symptoms include general fatigue, shortness of breath, swelling in the legs and ankles, and difficulty in accomplishing mild everyday activities such as walking or climbing stairs.^[3] Less than 50% of patients survive five years after initial diagnosis, and less than 25% of patients survive after 10 years of initial diagnosis.^[4] Though there is no cure, early diagnosis and treatment can allow patients to continue to maintain their quality of life and achieve a more normal life expectancy.

B-type natriuretic peptide (BNP) is a neurohormone secreted by the heart to regulate blood pressure and fluid balance.^[6] BNP levels in the blood plasma show strong correlation with the severity and prognosis of heart failure according to the New York Heart Association (NYHA) classification.^[6-8]

As CHF advances in severity, the level of BNP in blood plasma continues to increase.^[9,10] Additionally, the European Society of Cardiology(ESC) utilizes BNP testing in a diagnostic guideline for patient evaluation of heart failure.^[11]

3. Principle of the assay

The FREND BNP test cartridge is a one-time-use rapid “sandwich” immunoassay utilizing fluorescent nanoparticle in microfluidic flow to capture and quantify BNP in EDTA plasma specimens. A 35 μ L drop of patient EDTA plasma sample is placed in the FREND BNP cartridge, where the sample interacts with a proprietary mix of dry-loaded reagents. One of these reagents includes antibody-conjugated fluorescent nanoparticles, forming immune complexes with BNP in the patient samples. Capillary action moves the sample to the detection region, where capture antibodies grab the BNP-nanoparticle complexes. The fluorescence intensities of completely formed complexes are measured and BNP concentration is calculated by the FREND™ System.

4. Material provided

Q'ty	Contents	Catalogue number
25	FREND BNP Cartridge(s)	FRBN025
30	Disposable pipette tip(s)	
1	FREND BNP Code chip	
1	FREND BNP Package insert	

5. Materials required but not provided

- The FREND™ System
- Micro-pipette or an equivalent pipette capable of delivering 35 and 70 μ L
- Personal protective equipment and biohazard waste disposal containers

6. Warning and Precautions

- For *in vitro* diagnostic use only
- FREND BNP cartridges are single use.
- Use FREND Pipette and compatible tips, do not mouth pipette.
- Handle specimens in accordance with the OSHA Standard on Blood borne Pathogens.^[12]
- Wash hands thoroughly after handling specimens.
- Wear protective gloves, clothing, and eyewear.
- Heterophilic Antibody Interference: Some individuals have antibodies to mouse protein (HAMA), goat protein (HAGA), or rabbit protein (HARA); interferences may occur in FREND BNP.^[13,14]
- Do not use FREND BNP beyond indicated expiration date.
- Dispose of used FREND BNP cartridges and pipette tips as bio-hazardous materials.
- Do not attempt to disassemble FREND BNP cartridges: Micro-channels contain dry-loaded reagents that maybe bio-hazardous and/or toxic.
- Do not use FREND BNP cartridge if pouch is damaged or seal is broken.

7. Storage and Stability

Store FREND BNP cartridges at 2~8°C (35~46°F).

Do not remove the cartridge from sealed pouch until ready for use.

When stored properly, cartridges are stable until the marked expiration date.

Materials	Catalogue number
Refrigerator temperature storage (2~8 °C)	
FREND BNP cartridges	FRBN 025
Room temperature storage (18~25 °C)	
Pipette tips	None

8. Specimen collection and handling

- Human EDTA plasma sample (including EDTA plasma collected in separator tubes) are suitable for use with FREND BNP cartridges. Samples collected as serum in collection tubes or plasma in anticoagulant tubes (such as citrate or heparin) have not been validated for use in FREND BNP cartridges.
- Follow instructions detailed in this package insert as well as the specimen collection tube manufacturer's instructions for specimen collection and preparation (including manufacturer's instructions for centrifugation time and speed.)
- Collect blood by venipuncture technique into an EDTA plasma collection tube.
- After allowing the specimen to sufficiently mix with anticoagulant at room temperature, the sample tube can be centrifuged for 10 minutes at 3,000 rpm.
- The plasma should be immediately used.
- For optimal results, avoid grossly hemolytic, lipidic, or turbid specimens. Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter.
- When pipetting into the FREND BNP cartridge sample inlet, ensure that bubbles in the sample are avoided. Bubbles may restrict flow and result in an incomplete or erroneous test result.

9. Procedure

Code chip installation

There is no need for calibration to be performed by the end user. All calibration statistics and information have been electronically stored on the FREND BNP Code chip included in each box of FREND BNP cartridges. The FREND BNP Code chip is specific for each manufactured lot of FREND BNP cartridges.

Please refer to the FREND™ System user manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions follow here:

- 1) Press the 'Setup' button on the 'Main' screen.
- 2) Insert the FREND BNP Code chip into the code chip slot, following the direction arrow, at the rear of the FREND™ System.
- 3) Press the 'Code chip' button on the 'Setup' screen.
- 4) The information from the FREND BNP Code chip is automatically installed and saved.

- 5) When the installation is completed, press the 'OK' button to return to the 'Setup' screen.
- 6) Press the 'Item' button on the 'Setup' screen.
- 7) Confirm the FREND BNP cartridge lot number and code chip installation date. Press the 'Home' button to return to the 'Main' screen.

Quality control

• FREND™ System QC cartridges

FREND™ QC Cartridge contains multiple controls to check optic part of the system. By testing QC Cartridge, part of analytical components of the system of (1) laser power, (2) alignment, and (3) mechanical integrity are confirmed. For each day of patient testing, perform QC Cartridge testing. Refer to the quality control procedures section in the User Manual of FREND™ System. In brief, perform QC cartridge testing for the following conditions

- 1) Upon initial setup of the system,
- 2) Each day of patient testing,
- 3) When the system has been transported or moved,
- 4) Whenever there is uncertainty about the performance of the system,
- 5) Whenever required by your laboratory's quality control requirements.

• Internal procedural controls

FREND BNP cartridge contains built-in control features. Fluorescence signal in the Reference zone of each cartridge shows: (1) that enough volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this Reference zone signal is missing or lower than threshold, the FREND™ System consider it as an incorrect or failed test, not producing a test result but an error message. In addition, with each cartridge run, the system monitors, in part, for (1) flow of samples, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

• External quality control testing

Commercially available controls from a variety of manufacturers are available that contain BNP as a measured analyte. It is recommended that a minimum of at least two (2) levels of controls be run at least once per month or once for each new lot, whichever comes earlier. However, Controls should be run with a minimum frequency, depending on number of tests run in the

laboratory. Each laboratory should establish its own criteria based on the following parameters:

- 1) Each new lot,
- 2) Each new shipment (even if from the same lot previously received),
- 3) Each new operator (an individual who has not run the tests for at least two weeks),
- 4) Monthly, as a continued check on storage conditions,
- 5) Whenever problems (storage, operator, or other) are identified,
- 6) Or other times as required by your laboratory's standard QC procedures.

Individual laboratory policy will dictate exactly which control materials and lot numbers should be run, the frequency with which controls are to be tested, criteria for acceptance of the results and required corrective action to be taken if results do not meet laboratory criteria. If any external quality control sample values are out of the acceptable range, it will be necessary to investigate the problem before reporting patient results to assure there is not an instrument or software malfunction. Do not assay patient samples on the FRENDSTM System using FRENDS BNP if quality control results do not give expected values. Refer to your laboratory policies on how to determine acceptability of external control materials results. Each laboratory operates under a different set of regulations. Every laboratory must follow the standardized procedures acceptable to the regulatory agencies to which the laboratory is responsible.

Specimen processing

• Preparation

Remove sufficient FRENDS BNP cartridges from the refrigerator to test the number of patient samples and required external quality control materials. Allow the tubes and the sealed pouches containing the cartridges to come to room temperature for 15-30 minutes prior to the start of the testing sequence.

If using refrigerated patient samples, remove those from the refrigerator and allow them to come to room temperature prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed naturally and then mixed gently but thoroughly prior to testing. Testing should not begin on previously frozen samples until they have reached room temperature.

There are no other reagents or sample preparations necessary.

- **Specimen dilution procedure**

A specimen with BNP concentration greater than the linearity limit of the assay of 2,500 pg/mL may be diluted with a plasma that has been previously measured on the FREND BNP and found to contain < 30.0 pg/mL. The recommended dilution for samples with an initial result of >2,500.0 pg/mL is 1:9 (i.e. 1 mL of specimen with 9 mL of BNP low plasma). Dilute the sample so that the diluted sample reads between 200 and 2,000 pg/mL. Dilutions must be made manually and the final result on the diluted sample calculated manually by multiplying the result obtained on the diluted sample by the dilution factor (i.e. 10).

- **Assay procedure**

- 1) Open the FREND BNP box and remove a sealed FREND BNP pouch.
- 2) Allow cartridges to come to room temperature for 15–30 minutes prior to use.
- 3) Write patient information on the designated area on the FREND BNP cartridge.
- 4) Using the FREND Pipette, drop 35 µL of patient sample onto the indicated sample inlet of the FREND BNP cartridge.
- 5) In the 'Main' screen of the FREND™ System, press the 'Test' button on the touchscreen.
- 6) The system moves to the Patient ID screen automatically.
- 7) Type the Patient ID and press the 'Enter' button to begin the test.
- 8) Insert the cartridge into the cartridge slot using the cartridge arrow as a guide.

Caution: Please check the direction of the cartridge before insertion and assure the insertion is complete. It is recommended to insert the cartridge after the sample injection after 30 seconds elapsed in less than 5 minutes to obtain the optimal result of test.

- 9) When the reaction in the cartridge is completed, the FREND™ System will automatically begin the reading process.
- 10) When the measurements are completed, the cartridge will automatically be expelled and the results displayed.

Caution: Do not disconnect power cord or shut off power from the FREND™ System while a cartridge is in the reading chamber. This may cause a system error.

- 11) If the FREND™ System is connected to the optional printer, press the 'Print' button and the results will be output on the printer paper.

- 12) For more detailed instructions, please refer to the FREND™ System User Manual.

10. Procedural notes

If a specimen BNP concentration is found to be greater than the linearity limit of the assay of 2,500 pg/mL and a definitive result is required, the specimen should be diluted with low concentration sample that has been previously measured on the FREND BNP and found to contain <30 pg/mL BNP and then re-assayed according to the 'Assay Procedure'. The recommended dilution for samples with an initial result of >2,500 pg/mL is 1:20 (e.g., 10 µL of high level specimen and 190 µL of low level specimen). It is desirable to dilute the sample so that the dilute sample read between 100–2,000 pg/mL. Once a BNP result is obtained within the assay linearity, the original concentration of the unknown can be easily manually calculated.

Below is an example using a final diluted sample result of 253 pg/mL, a 1:20 dilution and a dilution sample concentration of 43 pg/mL where "X" equals the original concentration of the unknown sample and * indicates multiplication.

$$253 \text{ pg/mL} = \frac{(1 * X \text{ pg/mL}) + (19 * 43 \text{ pg/mL})}{20}$$

$$253 \text{ pg/mL} = \frac{X \text{ pg/mL} + 817 \text{ pg/mL}}{20}$$

$$(20) * (253) \text{ pg/mL} = X \text{ pg/mL} + 817 \text{ pg/mL}$$

$$5,060 \text{ pg/mL} = X \text{ pg/mL} + 817 \text{ pg/mL}$$

$$5,060 \text{ pg/mL} - 817 \text{ pg/mL} = X \text{ pg/mL}$$

$$4,243 \text{ pg/mL} = X \text{ pg/mL}$$

※ Original unknown concentration of BNP equals 4,243 pg/mL.

To use this formula to calculate the original concentration of any unknown sample BNP which exceeds the assay linearity and is diluted with another sample previously measured with the FREND BNP, just substitute the appropriate measured result on the diluted sample, the dilution factor and the diluent BNP concentration for those used above in the example. To do a quick check on your result, the original concentration for the unknown should always be just slightly less than the diluted result multiplied by the dilution factor without taking into consideration the BNP value of the sample used as the diluent.

11. Calculation of results

Results from FREND BNP are automatically displayed with Date/Time, User ID, Patient ID and BNP concentration on the screen of the FREND™ System.

The results are automatically saved and can be printed (with optional printer). For detailed information, refer to the FREND™ System User Manual.

Screen displayed for various concentration scenarios

Displayed result	Description
	BNP concentration Less than 30.0 pg/mL
	BNP concentration Not less than 30.0 pg/mL and not higher than 2,500.0 pg/mL
	BNP concentration Higher than 2,500.0 pg/mL

12. Limitations of the procedure

- 1) The results from FRENDBNP cartridge should be used as a diagnostic aid to the clinician, in conjunction with the patient's medical history, other laboratory tests, and clinical examination.
- 2) Specimens can be measured with the reportable range of 30.0 pg/mL ~ 2,500.0 pg/mL. Samples below the reportable range will yield "<30.0 pg/mL" and samples above the reportable range will yield ">2,500.0 pg/mL".
- 3) Specimens from patients with heterophilic antibodies, such as anti-mouse (HAMA), anti-goat (HAGA), or anti-rabbit (HARA) antibodies, maybe show falsely elevated or depressed values or may result in an incomplete test. Patients routinely exposed to animals or animal serum products can be prone to these types of heterophilic interferences.
- 4) FRENDBNP results should not be used interchangeably with other manufacturer's BNP or NT-proBNP test results.
- 5) FRENDBNP measurement should occur immediately before and 2 hours after the administration of recombinant BNP treatment (such as Nesiritide / Natrecor).

13. Performance characteristics

Precision

1) Repeatability

A precision study for FRENDBNP was performed as outlined in CLSI EP5-A3. Seven level standard samples were assayed in duplicates, twice per day over a period of 20 days.

Sample ID	Mean Conc. (pg/mL)	Within-run		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
A	0.00	3.115	23.1	1.630	12.1	0.702	5.2	3.586	26.6
B	50.00	3.970	7.9	1.271	2.5	1.051	2.1	4.299	8.6
C	100.00	6.316	6.4	1.319	1.3	1.301	1.3	6.582	6.7
D	200.00	9.864	4.9	1.240	0.6	16.227	8.1	19.030	9.5
E	400.00	21.697	5.4	5.979	1.5	4.202	1.1	22.895	5.7
F	1,000.00	54.426	5.5	17.689	1.8	7.120	0.7	57.669	5.8
G	2,500.00	130.873	5.2	48.465	1.9	22.043	0.9	141.288	5.7

2) Reproducibility (Lot)

Sample ID	Conc. (pg/mL)	Mean concentration (pg/mL)			Mean	SD	CV(%)
		197001	197002	197003			
A	0.0	13.99	10.57	13.16	12.57	1.79	14.20
B	50.0	47.72	53.45	50.92	50.70	2.87	5.66
C	100.0	95.56	108.19	97.64	100.46	6.77	6.74
D	200.0	186.00	200.65	189.68	192.11	7.62	3.97
E	400.0	405.30	401.51	395.49	400.77	4.95	1.23
F	1,000.0	945.81	950.01	1,000.20	965.34	30.27	3.14
G	2,500.0	2,508.48	2,565.30	2,519.81	2,531.20	30.07	1.19

3) Reproducibility (Operator)

Sample ID	Conc. (pg/mL)	Mean concentration (pg/mL)			Mean	SD	CV(%)
		Operator1	Operator2	Operator3			
A	0.0	12.22	11.31	11.08	11.54	0.60	5.20
B	50.0	47.64	50.85	49.35	49.28	1.60	3.25
C	100.0	92.74	105.12	101.75	99.87	6.40	6.41
D	200.0	205.85	200.66	203.04	203.18	2.60	1.28
E	400.0	395.97	404.78	427.05	409.27	16.01	3.91
F	1,000.0	1,006.88	984.26	1,060.64	1,107.26	39.23	3.86
G	2,500.0	2,479.84	2,483.93	2,513.66	2,492.48	18.46	0.74

Dilution linearity

To demonstrate the linearity of the assay, a serum base pool with an elevated BNP (4,000.0 pg/mL) was prepared and diluted to a total of 7 levels according to the dilution protocol outlined in CLSI EP6-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach. At each dilution level, the samples were tested in duplicate to determine the experimental value of BNP. Linearity was demonstrated from <30.0 pg/mL to >2,500.0 pg/mL. The measuring range for the FREND BNP is 30.0–2,500.0 pg/mL.


Analytical sensitivity

The limit of blank (LoB) and limit of detection (LoD) were determined according to the CLSI guideline EP17-A2. LoB was determined by 60 replicates of a BNP free sample. LoD was determined by 12 replicates of five low level patient samples.

LoB (Limit of Blank)	LoD (Limit of Detection)	LoQ (Limit of Quantitation)
20.42 pg/mL	26.80 pg/mL	26.80 pg/mL

Comparative analysis

FREND BNP was compared to the Predicate device using guidelines outlined in CLSI document EP9-A2. Clinical samples (EDTA plasma, n=120) were measured in duplicate on both systems. Linear regression analysis demonstrated a slope of 1.0338 and a correlation coefficient (r) of 0.983.

Cross reactivity

The following substances were evaluated for potential cross-reactivity with the FREND BNP at three concentrations. Testing was done according to the CLSI protocol EP7-A2. No significant cross-reactivity was found.

Substance	Concentration (pg/mL)
ANP	1,000
CNP	1,000
NT-proBNP	1,000
Angiotensin I	600
Angiotensin II	600
Angiotensin III	1,000

Interference

1) Endogenous

The interference study was performed as recommended in the CLSI EP7-A2 protocol using three concentrations of BNP. Recovery within 90% to 110% of the expected BNP was considered as lack of interference. No interference by the substances below was found.

No.	Substances	Concentration
1	Hemoglobin	500 mg/dL
2	Bilirubin	20 mg/dL
3	Triglycerides	3 g/dL
4	Albumin	3 g/dL

2) Drug

A variety of common medications were tested for interference with the FRENED BNP. The testing showed there was no significant interference (<10%) from the tested drugs that would affect the interpretation of a Brain Natriuretic Peptide result as assayed on the FRENED BNP.

No.	Substances	Concentration
1	Acetaminophen	30 µg/mL
2	Atenolol	10 µg/mL
3	Caffeine	60 µg/mL
4	Captopril	5 µg/mL
5	Chloramphenicol	50 µg/mL
6	Diclofenac	50 µg/mL
7	Digoxin	2 ng/mL
8	Diltiazem	40 µg/mL
9	Enalapril Maleate	300 µg/mL
10	Hydrochlorothiazide	6 µg/mL
11	Indomethacin	36 µg/mL
12	Nifedipine	400 µg/mL
13	Probenecid	600 µg/mL
14	Propranolol	2 µg/mL
15	Quinidine	12 µg/mL
16	Sulfamethoxazole	400 µg/mL
17	Trimethoprim	40 µg/mL

14. References


- 1) American Heart Association, Heart Disease and Stroke Statistics – 2012 update. *Circulation*. 2012, 125: e2-e220.
- 2) http://www.cdc.gov/dhdsdp/data_statistics/fact_sheets/fs_heart_failure.htm
- 3) Wu, A. H., B-type Natriuretic Peptide and its Clinical Utility in Patients with Heart Failure. *MLO Medical Laboratory Observer*. 2001, 10: 10-14.
- 4) Loh, J.C. et al. Temporal Trends in Treatment and Outcomes for Advanced Heart Failure with Reduced Ejection Fraction from 1993-2010; Findings from a University Referral Centre. *Circulation Heart Fail*. 2013, 6: 411-419.
- 5) Mukoyama, M. et al. Brain natriuretic peptide as a novel cardiac hormone in humans: Evidence for an exquisite dual natriuretic peptide system, atrial natriuretic peptide and brain natriuretic peptide. *Journal of Clinical Investigation*. 1991, 87: 1402-1412.
- 6) Maisel, A.S. et al. Utility of B-natriuretic peptide as a rapid, point-of-care test for screening patients undergoing echocardiography to determine left ventricular dysfunction. *American Heart Journal*. 2001, 141: 367-374.
- 7) Wallen, T. et al. Brain natriuretic peptide predicts mortality in the elderly. *Heart*. 1997, 77: 264-267.
- 8) Wiecezorek, S.J. et al. A rapid B-type natriuretic peptide assay accurately diagnoses left ventricular dysfunction and heart failure a multicenter evaluation. *American Heart Journal*. 2002, 144: 834-839.
- 9) Yasue, H. et al. Localization and mechanism of secretion of B-type natriuretic peptide in comparison with those of A-type natriuretic peptide in normal subjects and patients with heart failure. *Circulation*. 1994, 90: 195-203.
- 10) Clerico, A. et al. Circulating levels of cardiac natriuretic peptides (ANP and BNP) measured by highly sensitive and specific immunoradiometric assay s in normal subjects and in patients with different degrees of heart failure. *Journal of Endocrinological Investigation*. 1998, 21: 170-179.
- 11) Schroff, R.W. et al. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45: 879-885.
- 12) Chao, E.L., Henshaw, J.L. Occupational Safety and Health Administration: Model Plans and Programs for the OSHA Bloodborne Pathogens and Hazard Communications Standards. OSHA 3186-06R, 2003.
- 13) Schroff, R.W., Foon, K.A., Beatty, S.M., Oldham, R.K., Morgan, A.C. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45: 879-885.
- 14) Boscatto, L.M. et al. Heterophilic Antibodies: A Problem for All Immunoassay. *Clinical Chemistry*. 1988, 34: 27.

Glossary of Symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
 www.nanoentek.com/eifu.php	Consult Instructions for Use An electronic instructions for use (eIFU) indicator (website address) may accompany the symbol when used to indicate an instruction to consult an eIFU.
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only CAUTION: Federal (U.S.) law restricts this device to sale by or on order of a physician.
	US Corporation
	Patient ID
	Result
 Sample Drop	Sample Drop
	Authorized representative in the European Community

UK Representative	Authorized representative in United Kingdom
CH REP	Authorized representative in Switzerland
BRH	Authorized representative in Brazil

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14. Riferimenti

- 1) American Heart Association, Heart Disease and Stroke Statistics – 2012 update. *Circulation*. 2012, 125: e2-e220.
- 2) http://www.cdc.gov/dhdsdp/data_statistics/fact_sheets/fs_heart_failure.htm
- 3) Wu, A. H., B-type Natriuretic Peptide and its Clinical Utility in Patients with Heart Failure. *MLO Medical Laboratory Observer*. 2001, 10: 10-14.
- 4) Loh, J.C. et al. Temporal Trends in Treatment and Outcomes for Advanced Heart Failure with Reduced Ejection Fraction from 1993-2010; Findings from a University Referral Centre. *Circulation Heart Fail*. 2013, 6: 411-419.
- 5) Mukoyama, M. et al. Brain natriuretic peptide as a novel cardiac hormone in humans: Evidence for an exquisite dual natriuretic peptide system, atrial natriuretic peptide and brain natriuretic peptide. *Journal of Clinical Investigation*. 1991, 87: 1402-1412.
- 6) Maisel, A.S. et al. Utility of B-natriuretic peptide as a rapid, point-of-care test for screening patients undergoing echocardiography to determine left ventricular dysfunction. *American Heart Journal*. 2001, 141: 367-374.
- 7) Wallen, T. et al. Brain natriuretic peptide predicts mortality in the elderly. *Heart*. 1997, 77: 264-267.
- 8) Wiecezorek, S.J. et al. A rapid B-type natriuretic peptide assay accurately diagnoses left ventricular dysfunction and heart failure a multicenter evaluation. *American Heart Journal*. 2002, 144: 834-839.
- 9) Yasue, H. et al. Localization and mechanism of secretion of B-type natriuretic peptide in comparison with those of A-type natriuretic peptide in normal subjects and patients with heart failure. *Circulation*. 1994, 90: 195-203.
- 10) Clerico, A. et al. Circulating levels of cardiac natriuretic peptides (ANP and BNP) measured by highly sensitive and specific immunoradiometric assay s in normal subjects and in patients with different degrees of heart failure. *Journal of Endocrinological Investigation*. 1998, 21: 170-179.
- 11) Schroff, R.W. et al. Human Anti-murine immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45: 879-885.
- 12) Chao, E.L., Henshaw, J.L. Occupational Safety and Health Administration: Model Plans and Programs for the OSHA Bloodborne Pathogens and Hazard Communications Standards. OSHA 3186-06R, 2003.
- 13) Schroff, R.W., Foon, K.A., Beatty, S.M., Oldham, R.K., Morgan, A.C. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45: 879-885.
- 14) Boscatto, L.M. et al. Heterophilic antibodies: A Problem for All Immunoassay. *Clinical Chemistry*. 1988, 34: 27.

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NanoEntek, Inc.
MT Promedt Consulting GmbH