


FDA 510K Cleared



Testosterone



Immediate Consultation Testosterone



FAST



EASY TO USE



ACCURATE



FRENDS™ System

Quantitative immunoassay analyzer
for Testosterone hormone

- Improve men's healthcare
- Increase customer satisfaction
- Increase practice efficiency
- Maximize your profit
- Test range 20 – 1,500 ng/dL

FREND™ System

Accelerating decision-making

Accurate and fully automated FREND™ System brings fast results. With a simple 3-step operation and a wide range of disease markers, the system supports quick decision-making.




Simple and Fast testing



Once a month liquid quality control



User-friendly Icon



6 seconds code chip calibration

Test procedures



STEP 01
Add sample



STEP 02
Insert AP tube



STEP 03
Get Result




Testosterone

Specifications

Assay method	Fluorescent immunoassay
Analytical measurement range	20 ng/dL ~ 1,500 ng/dL
Sample type	Serum or Plasma (Lithium heparin or EDTA)
Sample volume	70 µL
Limit of detection	20 ng/dL
Reaction time	Cartridge reaction: < 3 min.
Pre-incubation	AP (5 minutes)
Package	20 tests/kit
Storage condition	2°C – 8 °C (35°F - 46°F)
Shelf life	18 months

Correlation



Ordering information

FREND™ System	Cat# F10
FREND™ Testosterone	Cat# FRTEAP 020
CPT Code	84403



NanoEntek [website](http://www.nanoentek.com) www.nanoentek.com
e-mail ivdst@nanoentek.com

NESCT-FRTEAP-002EN(V.0.1)

NanoEntek, Inc.

Head Office
12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea
Tel +82-2-6220-7940 / Fax +82-2-6220-7999

NanoEntek America, Inc.
220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA
Tel +1-781-472-2558 / Tel CS: +1-781-216-8916 / Fax: +1-781-790-5649

NanoEntek Europe | med-tech supplies GmbH
Lochhamerstr. 4a, 82152 Martinsried, Germany
Tel +49-89-21-55-38-43 / Fax +49-89-99-95-46-60